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Bibliografia

- [1] Holden BA, Fricke TR, Wilson DA, Jong M, Naidoo KS, Sankaridurg P, *et al.* Global Prevalence of Myopia and High Myopia and Temporal Trends from 2000 through 2050. *Ophthalmology*. 2016;123:1036-42
- [2] McCullough SJ, O'Donoghue L, Saunders KJ. Six Year Refractive Change among White Children and Young Adults: Evidence for Significant Increase in Myopia among White UK Children. *PLoS One*. 2016;11:e0146332
- [3] Dong L, Fazzari M, Gwiazda J, Comet. Myopia stabilization and associated factors among participants in the Correction of Myopia Evaluation Trial (COMET). *Invest Ophthalmol Vis Sci*. 2013;54:7871-84
- [4] Donovan L, Sankaridurg P, Ho A, Naduvilath T, Smith EL, 3rd, Holden BA. Myopia progression rates in urban children wearing single-vision spectacles. *Optom Vis Sci*. 2012;89:27-32
- [5] Wolffsohn JS, Calossi A, Cho P, Gifford K, Jones L, Li M, *et al.* Global trends in myopia management attitudes and strategies in clinical practice. *Cont Lens Anterior Eye*. 2016;39:106-16
- [6] Cho P, Cheung SW. Retardation of Myopia in Orthokeratology (Romio) Study: A 2 year Randomized Clinical Trial. *Invest Ophthalmol Vis Sci*. 2012;53:7077–85.
- [7]. Cho P, Cheung SW, Edwards M. The Longitudinal Orthokeratology Research in Children (Loric) in Hong Kong: A Pilot Study on Refractive Changes and Myopic Control. *Curr Eye Res* 2005;30:71–80
- [8]. Walline JJ, Jones LA, Sinnott LT. Corneal Reshaping and Myopia Progression. *Br J Ophthalmol* 2009;93:1181–5
- [9]. Anstice NS, Phillips JR. Effect of Dual-focus Soft Contact Lens Wear on Axial Myopia Progression in Children. *Ophthalmology* 2011;118:1152–61
- [10]. Sankaridurg P, Holden B, Smith E, 3rd, *et al.* Decrease in Rate of Myopia Progression with a Contact Lens Designed to Reduce Relative Peripheral Hyperopia: One-year Results. *Invest Ophthalmol Vis Sci*. 2011; 52:9362–7
- [11]. Walline JJ, Greiner KL, McVey ME, *et al.* Multifocal Contact Lens Myopia Control. *Optom Vis Sci*. 2013;90: 1207–14
- [12]. Fujikado T, Ninomiya S, Kobayashi T, *et al.* Effect of Low-addition Soft Contact Lenses with Decentered Optical Design on Myopia Progression in Children: A Pilot Study. *Clin Ophthalmol*. 2014;8:1947–56
- [13]. Lam CS, TangWC, Tse DY, *et al.* Defocus Incorporated Soft Contact (Disc) Lens Slows Myopia Progression in Hong Kong Chinese Schoolchildren: A 2-year Randomised Clinical Trial. *Br J Ophthalmol* 2014;98:40–5
- [14]. Paune J, Morales H, Armengol J, *et al.* Myopia Control with a Novel Peripheral Gradient Soft Lens and Orthokeratology: A 2-year Clinical Trial. *Biomed Res Int* 2015;2015:507572

- [15]. Aller TA, Liu M, Wildsoet CF. Myopia Control with Bifocal Contact Lenses: A Randomized Clinical Trial. *Optom Vis Sci* 2016;93:344–52
- [16]. Cheng X, Xu J, Chehab K, *et al.* Soft Contact Lenses with Positive Spherical Aberration for Myopia Control. *Optom Vis Sci* 2016;93:353–66
- [17]. Ruiz-Pomeda A, Perez-Sanchez B, Valls I, *et al.* MiSight Assessment Study Spain (MASS). A 2-Year Randomized Clinical Trial. *Graefes Arch Clin Exp Ophthalmol* 2018;256:1011–21
- [18]. Allen PM, Radhakrishnan H, Price H, *et al.* A Randomised Clinical Trial to Assess the Effect of a Dual Treatment on Myopia Progression: The Cambridge Anti-Myopia Study. *Ophthalmic Physiol Opt* 2013; 33:267–76
- [20] Smith EL, 3rd. Prentice Award Lecture 2010: A case for peripheral optical treatment strategies for myopia. *Optom Vis Sci.* 2011;88:1029-44
- [21] Smith EL, 3rd, Hung LF, Arumugam B. Visual Regulation of Refractive Development: Insights from Animal Studies. *Eye (Lond)* 2014;28:180–8
- [22] Liu Y, Wildsoet C. The Effect of Two-zone Concentric Bifocal Spectacle Lenses on Refractive Error Development and Eye Growth in Young Chicks. *Invest Ophthalmol Vis Sci* 2011;52:1078–86
- [23] Liu Y, Wildsoet C. The Effective Add Inherent in 2-zone Negative Lenses Inhibits Eye Growth in Myopic Young Chicks. *Invest Ophthalmol Vis Sci.* 2012;53: 5085–93
- [24] Arumugam B, Hung LF, To CH, *et al.* The Effects of Simultaneous Dual Focus Lenses on Refractive Development in Infant Monkeys. *Invest Ophthalmol Vis Sci* 2014;55:7423–32
- [25] Benavente-Perez A, Nour A, Troilo D. The Effect of Simultaneous Negative and Positive Defocus on Eye Growth and Development of Refractive State in Marmosets. *Invest Ophthalmol Vis Sci* 2012;53:6479–87
- [26] Arumugam B, Hung LF, To CH, *et al.* The Effects of the Relative Strength of Simultaneous Competing Defocus Signals on Emmetropization in Infant Rhesus Monkeys. *Invest Ophthalmol Vis Sci* 2016;57: 3949–60
- [27] Tse DY, Lam CS, Guggenheim JA, *et al.* Simultaneous Defocus Integration During Refractive Development. *Invest Ophthalmol Vis Sci* 2007;48: 5352–9
- [28] Chamberlain P. *et al.* A 3-Year Randomized Clinical Trial of MiSight Lenses for Myopia Control. *Optom Vis Sci.* 2019;96:556–567 doi:10.1097/OPX.0000000000001410
- [29] Flitcroft DI. The complex interactions of retinal, optical and environmental factors in myopia aetiology. *Prog Retin Eye Res.* 2012;31:622-60
- [30] Rose K, Harper R, Tromans C, Waterman C, Goldberg D, Haggerty C, *et al.* Quality of life in myopia. *Br J Ophthalmol.* 2000;84:1031-4
- [31] Dias L, Manny RE, Weissberg E, Fern KD. Myopia, contact lens use and self-esteem. *Ophthalmic Physiol Opt.* 2013;33:573-80
- [32] Walline JJ, Gaume A, Jones LA, Rah MJ, Manny RE, Berntsen DA, *et al.* Benefits of contact lens wear for children and teens. *Eye Contact Lens.* 2007;33:317-21
- [33] Walline JJ, Jones LA, Sinnott L, Chitkara M, Coffey B, Jackson JM, *et al.* Randomized trial of the effect of contact lens wear on self-perception in children. *Optom Vis Sci.* 2009;86:222-32
- [34] Walline JJ, Long S, Zadnik K. Daily disposable contact lens wear in myopic children. *Optom Vis Sci.* 2004;81:255-9
- [35] Walline JJ, Lorenz KO, Nichols JJ. Long-term contact lens wear of children and teens. *Eye Contact Lens.* 2013;39:283-9
- [36] Bullimore MA. The Safety of Soft Contact Lenses in Children. *Optom Vis Sci.* 2017;94:638-46
- [37] Chalmers RL, Wagner H, Mitchell GL, Lam DY, Kinoshita BT, Jansen ME, *et al.* Age and other risk factors for corneal infiltrative and inflammatory events in young soft

contact lens wearers from the Contact Lens Assessment in Youth (CLAY) study. Invest Ophthalmol Vis Sci. 2011;52:6690-6

[38] Chamberlain P. *et al.* Further comparison of myopia progression in new and established myopia control treatment (MiSight® groups). BCLA Clinical Conference 2019

[39] Cho P, Boost M. Daily disposable lenses: the better alternative. Contact Lens Anterior Eye. 2013;36:4-12

[40] Chalmers RL, Hickson-Curran SB, Keay L, Gleason WJ, Albright R. Rates of adverse events with hydrogel and silicone hydrogel daily disposable lenses in a large postmarket surveillance registry: the TEMPO Registry. Invest Ophthalmol Vis Sci. 2015;56:654-63